
# ENA SHINY

“Natural Technology”

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CE 0123

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## (EN) FINISHING KIT, ENA SHINY

A modern esthetic composite restoration foresees the accurate preparation of the cavity, the use of a stratification technique through dentin fluorescent masses and opalescent enamel masses and an accurate finishing and polishing phase, able to give the restoration surface a similar feature to that of a natural tooth.

### **Rotative instruments**

The rotative instruments used in the finishing and polishing of composite restorations consist of diamond burs and rubbers. Diamond burs have different grain measures and are usually used in sequences: those with a bigger grain are used for the carving and first working of the restoration, those with a smaller grain for the finishing. Rubbers can be used exclusively in the last phase of finishing, being aware not to push too much on the tooth surface, in order not to overheat it, creating pulpar damage, and in order not to burn the resin matrix of the composite material.

### **Abrasive discs and strips**

The clinic use of flexible discs is limited because they are not able to reach certain anatomic areas, such as ditches, and they generally reduce the surface characterisations, hardly achieved during the carving phase of restoration. Therefore their use must be limited to incisal and interproximal areas. Abrasive strips are used for the finishing and polishing of interproximal edges. Great care must be taken during their use, in order not to damage the parodontal legament and the gingival papilla and in order to respect the holding of an interproximal contact point. These are used in sequences as well and foresee different passages to obtain a good final result.

### **Polishing pastes**

Hybrid composites of the last generation respond to polishing very well, through diamond pastes and aluminium oxide pastes. The Shiny pastes are available in two different grain measures: 3 microns (paste A) and 1 micron (paste B), to be used in sequence with brushes. A final polishing is achieved with felts, using the aluminium oxide paste (paste C).

### **CLINICAL PROCEDURE**

**Pic. 1-2.** Before starting the finishing of a composite reconstruction, it is necessary to apply an Air Block (Shiny G), when the restoration is finished and before the final light curing. This glycerin-based product eliminates the oxygen from the surface (that inhibits the complete curing of the composite) and increases the resistance of the surface material. Tests have shown that using an air block in combination with diamond pastes the surface hardness increases of 10% and its deviation from medium value decreases from 1,2% to 0,8%, thus the surface is more even. (**Pic. 2** dev=hardness deviation from medium value, AB= Air Block, DP= polishing with diamond paste).

**Pic. 3-4.** The anatomy of an anterior tooth is often characterised by a rich surface morphology, represented by ditches and lumps (macrogeography), and by growth lines of the enamel (microgeography). The first finishing phase of a composite restoration foresees the use of Ena Shiny diamond burs in a fine grain and/or multiblade burs. With these burs, used at a low speed, a revision of the surface morphology, obtained through carving, is made; the final macrogeography is a result. In order to facilitate this phase, the obtained surface can be highlighted, during carving and finishing (macrogeography), by using the Shiny L varnish. The color subtraction technique allows to clearly visualize the depth of ditches and the disposition of lumps.

**Pic. 5.** In interproximal areas, the work started with the hand instruments can be beautifully finished with diamond strips (available in fine grain red and extra-fine yellow, in 2 and 4 mm height) or

1. Air block applied during the final polymerisation

### SURFACE TREATMENT

	VICKERS H.	DEV.
ENA Composite	55	1,2
ENA Composite (AB + DP)	62	0,8
MICROFILL	18-30	4


2. Composite surface treatment: Vickers hardness and its deviation (AB=air block, DP=diamond paste)

3. Macrogeography can be corrected by using fine grain diamond burs (red ring: 861514012; 853514008)...

4. ... or a multi-blade bur (CS9904).


5. Diamond and paper strips and mini-discs for interproximal area

6. Diamond silicon rubber (Shiny 14) for pre-polishing

7. Microgeography can be checked by marking with pencil ....

8. ... and is enhanced with a carbundum bur (Shiny 21) or a multi-blade (CS9904)

9. A 3 micron diamond paste (Shiny A) is applied...

10. ... with a soft hair goat brush (Shiny S)


11. A 1 micron diamond paste (Shiny B) is applied...


12. ... with a soft hair goat brush (Shiny S)

13. Hi-polishing with Al-Oxide (Shiny C)


14. ... and a felt (Shiny F or SF)

15. Only convex areas (marked with stars) need to be polished with aluminium oxide paste (Shiny C)

16. Final polished composite restoration

17. Impregnated brushes (silicon carbide)/ hair goat brush mini

#### EKSKIT SHINY FINISHING KIT

- SHINY A diamond paste 3 micron 2 g
- SHINY B diamond paste 1 micron 2 g
- SHINY C shiny paste 35 g • SHINY G air block 7 tips
- SHINY L varnish 3 ml • SHINY CS burs 4 pcs
- SHINY 14 rubber point 2 pcs • SHINY S/SS hair goat brush 3 pcs
- SHINY F/SF felt wheel 5 pcs • SHINY FD felt disc 5 pcs
- SHINY D... shiny disc 12 pcs • SHINY 5...shiny strip 4 pcs
- SHINY M... shiny strip 12 pcs • Mandrel 2 pcs
- SHINY 4... brush for posterior 3 pcs
- SHINY 1 Handle with 10 brushes (SHINY 50)

#### ES+G POLISHING PASTE SHINY KIT

- SHINY A diamond paste 3 micron 2 g
- SHINY B diamond paste 1 micron 2 g
- SHINY C shiny paste 35 g • SHINY G air Block 2 tips
- SHINY 14 rubber point 1 pc • SHINY S/SS hair goat brush 3 pcs
- SHINY F/SF felt wheel 3 pcs • Mandrel 1 pc
- SHINY FD felt disc 3 pcs
- SHINY 5... shiny strip 4 pcs
- SHINY M... shiny strip 4 pcs
- SHINY 4... brush for posterior 2 pcs

For contents of other Shiny Kits see on the package.

18. Ena Shiny Kit

paper abrasive strips (Shiny SM and SMS, with double grains, medium-coarse, with two widths, Shiny SF and SFS with double grains fine and super-fine, with two widths) and discs (Shiny DM medium and DF super-fine) reducing the build up of bacterial plaque in these areas that are difficult to access.

**Pic. 6.** Finishing is accomplished by using a diamond silicone rubber that has the shape of a flame (Shiny 14), dunked in pumice. This instrument, used at a low speed and under a spray of water, levigates the surface without erasing the macrogeography.

**Pic. 7-8** The microgeography is created with the help of a Shiny carborundum bur (Shiny 21) o multi-blade bur (CS9904), which should be used with great care. By rubbing the dental surface with an articulation paper, the intensity of the created macrogeography is pointed out.

**Pic. 9-10-11-12-17.** The pre-finishing, obtained with the silicone rubber, prepares the restoration to a proper finishing, that makes the surface shiny, without destroying the micro and macrogeography. For this, a goat brush (Shiny S for anterior, Shiny SS for posterior) is used, together with diamond pastes Shiny A and B, with different grains (3 microns and 1 micron). The surface must first be treated without water and at a very low speed, with short working times, following the inclination illustrated in the picture, in order to avoid over-heat; afterwards, a heavy spray of water may be used, increasing speed slightly. For posterior, Shiny 41-42-43 brushes impregnated with silicone carbide can be used (Pic.17)

**Pic. 13-14-15.** Convex areas of the surface, more exposed to using, always appear shinier than the concave ones; therefore with a round felt (Shiny F or Shiny SF) or a disc felt (Shiny FD) and a paste Shiny C in aluminium oxide, we obtain a brighter reflected light in these areas. The paste can be used to polish composite restorations every time a patient returns for prophylaxis (maintaining phase).

**Pic. 16.** If made properly, these different working phases that can be done in one sitting, allow the clinician to obtain a restoration that really simulates the real tooth, with a guarantee of the respect of mechanical characteristics of used material, dental tissues and parodontal tissues.

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**Side effects and contra-indications**

- Do not use discs on the vestibular surfaces in order to avoid flattening the surface (they can be used interproximally and on the incisal edges).
- Use rotating instruments at low speed to avoid overheating the pulp (or the tooth structure - risk of necrosis) and the composite (risk to burn the resin matrix).

**Take care:**

- SHINY L: Shake well before use. Close cap tightly after use.
- A-B-C-G PASTE: If the material has not been used for a long time, dispose of the dried-out material at the tip of the syringe. Close cap tightly after use. MSDS available on web site [www.enahri.com](http://www.enahri.com)

**Note:** instruments and materials do not contain hazardous components. They are safe to use when instructions are followed.

**Precautions:**

- Use rotating instruments at low speed and intermittently. Use in conjunction with water spray to reduce heat production.
- Diamond Strips, SHINY41-42-43 impregnated brushes for posterior, diamond silicone rubber are reusable. SHINY S and SHINY SS brushes and other instruments are disposable.
- Use spray with suction to remove glycerin and varnish.

## REPROCESSING INSTRUCTIONS FOR REUSABLE MEDICAL DEVICES (ISO EN 17664)

### WARNINGS

Instruments are supplied not sterile and should be sterilized before the first use and after each use. Appropriate handling precautions should be taken during use and reprocessing. Gloves, eye protection and a mask should be worn.

### Limitations on reprocessing

Repeated processing has minimal effect on these instruments; the working part is deteriorating mainly with use.

### INSTRUCTIONS

**Point of use:** Remove excess soil under running water.

**Container:** No particular requirements. They should be protected from damage to the working part: burs holder is recommended.

**Transportation:** It is recommended that instruments are reprocessed as soon as possible after using them. If transported wet there is an increased chance of staining or corrosion. Prolonged storage in disinfectant solutions may result in corrosion and should be avoided.

**Preparation for cleaning:** Remove the screw and the rubber from the head of the mandrel eliminating the abrasive wheel mounted (only for mandrel).

### Cleaning:

- **automated cleaning:** Ultrasonic Cleaner as Biosonic with general purpose liquid as UC30 or other liquids validated for cleaning of dental rotary instruments (manufacturers' instructions must be followed). Specific disinfectant solutions for ultrasonic can also be used only if validated for use with rotating instruments, specifically, for brushes and rubber points. 1. Load burs directly in a beaker or in a burs holder. 2. Run cycle, minimum 15 minutes. 3. When unloading check instruments and tips for complete removal of visible soil. If necessary repeat cycle or use manual cleaning.

- **manual cleaning:** 1. Rinse excess soil from instrument under running cold water. 2. keep the instruments immersed in a sink reserved for the purpose in which case a disinfectant agent validated for processing of dental rotary instruments must be used and the disinfectant manufacturers' instructions must be followed and brush and use wire brushes to remove any contamination Care should be taken to avoid spreading contaminants by spraying or splashing during the brushing process 3. Rinse under clean running water and than better again with ionized water.

**Disinfection:** Disinfectant solution may be used in accordance with label instructions. Only disinfectant agent validated for processing of dental rotary instruments must be used.

**Drying:** When drying is achieved as part of a washer disinfector cycle do not exceed 140°C.

**Maintenance:** No particular requirements.

**Inspection and Function - Testing:** After clearing, inspect all instruments, with the aid of magnification if necessary, to ensure that all contamination has been removed. Repeat the cleaning process if necessary. Visually inspect for damage and eliminate broken or deteriorated instruments.

**Packaging:** Single: A standard packaging material may be used. In sets: Instruments may be loaded into dedicated bur holder, Wrap bur holder using appropriate method.

**Sterilization:** Vacuum autoclave (ISO 17665). Cycle: 3 minutes minimum (18 minutes for diamond and carbide burs and Ena Post drills), 134°C maximum, 2,1 bar (follow autoclave manufacturers instructions).

**Storage:** No specific limitations related to storing duration. Storage should be in dry, clean conditions and at ambient temperature. Containers or pouches must be dry before opening to avoid recontamination of the contents with water.

**Additional Information:** When sterilizing multiple instruments in one autoclave cycle ensure that the sterilizer's maximum load, suggested by the manufacturer, is not exceeded.

**Manufacturer contact:** See brochure for telephone and address of local representative or telephone (+39) 0185 7887880

The instructions provided above have been validated by the medical device manufacturer as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using equipment, materials and personnel in the processing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

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2. Tratamiento de superficie: dureza, Vickers y su desviación (AB= air block, DP= pasta diamantada