

ENA MATRIX

Sectional Matrix "THIN" 35 micron
Ideal when there are space problems



EMKC



EMC35



EMC45



EMC55



EMC65

ENA WEDGE

Available 3 sizes of plastic wedges and 4 sizes of wooden wedges

- They can anatomically adapt to the interproximal area and allow a complete seal both on the vestibular and lingual side
- In plastic wedges a hole, in the handle, allows the positioning and removal with the specific tweezers without any problem.



EMWW



EMWWXS



EMWWS



EMWWM



EMWWL

Instruments



EMAS

Ref	Description	Q.ty
EMKC	SECTIONAL MATRIX KIT "THIN" 35 MICRON	
	Sectional matrix 35 micron Small	25 pcs
	Sectional matrix 35 micron Medium	25 pcs
	Sectional matrix 35 micron Large	25 pcs
	Sectional matrix 35 micron Subgingival	25 pcs
	Extra small wooden wedge (orange)	25 pcs
	Small wooden wedge (white)	25 pcs
	Medium wooden wedge (yellow)	25 pcs
	Large wooden wedge (pink)	25 pcs
EMC35	Sectional matrix 35 micron Small	50 pcs
EMC45	Sectional matrix 35 micron Medium	50 pcs
EMC55	Sectional matrix 35 micron Large	50 pcs
EMC65	Sectional matrix 35 micron Subgingival	50 pcs

Plastic wedge

Ref	Description	Q.ty
EMW0S	Small plastic wedge (white)	100 pcs
EMW0M	Medium plastic wedge (pink)	100 pcs
EMW0L	Large plastic wedge (violet)	100 pcs

Wooden wedge

Ref	Description	Q.ty
EMWW	WOODEN WEDGE KIT	
	Extra small wooden wedge (orange)	100 pcs
	Small wooden wedge (white)	100 pcs
	Medium wooden wedge (yellow)	100 pcs
	Large wooden wedge (pink)	100 pcs
EMWWXS	Extra small wooden wedge (orange)	100 pcs
EMWWS	Small wooden wedge (white)	100 pcs
EMWWM	Medium wooden wedge (yellow)	100 pcs
EMWWL	Large wooden wedge (pink)	100 pcs

Ref	Description	Q.ty
EMAS	Stainless steel ring with "V" shape tines	1 pc