

Microhybrid, radiopaque light curing composite

"HFO" the classical composite for direct and indirect aesthetic reconstructions.

- Fluorescent dentine calibrated to the natural tooth
- Generic Enamels: three values to reproduce elderly (GE1), adult (GE2) and young (GE3) teeth
- Opalescent and Intensive enamels for reproducing opalescent and characterization of the natural enamel

Ref	Description	Syringes
CHFGE1	Enamel plus HFO generic enamel GE1	5 g
CHFGE2	Enamel plus HFO generic enamel GE2	5 g
CHFGE3	Enamel plus HFO generic enamel GE3	5 g
CHFOW	Enamel plus HFO opalescent enamel White	5 g
CHFOG	Enamel plus HFO opalescent enamel Grey	5 g
CHFOA	Enamel plus HFO opalescent enamel Amber	5 g
CHFOBN	Enamel plus HFO opalescent Blue natural	5 g
CHFIM	Enamel plus HFO intensive Milky	5 g
CHFIW	Enamel plus HFO intensive White	5 g
CHFUD1	Enamel plus HFO dentine UD1 (A1)	5 g
CHFUD2	Enamel plus HFO dentine UD2 (A2)	5 g
CHFUD3	Enamel plus HFO dentine UD3 (A3)	5 g
CHFUD3,5	Enamel plus HFO dentine UD3,5 (A3,5)	5 g
CHFUD4	Enamel plus HFO dentine UD4 (A4)	5 g
CHFUD5	Enamel plus HFO dentine UD5	5 g
CHFUD6	Enamel plus HFO dentine UD6	5 g
COSHFOB1	Enamel plus HFO dentine B1	5 g
COSHFOB2	Enamel plus HFO dentine B2	5 g
COSHFOB3	Enamel plus HFO dentine B3	5 g
COSHFOB4	Enamel plus HFO dentine B4	5 g
COSHFOC1	Enamel plus HFO dentine C1	5 g
COSHFOC2	Enamel plus HFO dentine C2	5 g
COSHFOC3	Enamel plus HFO dentine C3	5 g
COSHFOC4	Enamel plus HFO dentine C4	5 g
COSHFOC3	Enamel plus HFO dentine D3	5 g

* Colours B...C...D3 shades available on request (colours of Vita® shade guide. Vita® is a registered trademark of Vita Zahnfabrik H. Rauter mbH & Co. KG, Bad Säckingen - D)

*Available also in unidose Tips Kit (Ref. CHF15T) and 12 pcs refill (Ref. CHFTGE1, CHFT...)



Ref	Description	Q.ty
CHF15	ENAMEL PLUS HFO Complete Kit including: 15 syringes 5 g with holder 7 dentine: UD1 (A1), UD2 (A2), UD3 (A3), UD3,5 (A3,5), UD4 (A4), UD5, UD6 3 generic enamel: GE1, GE2, GE3 3 opalescent enamel: OBN, OA, OW 2 intensive enamel: White (IW), Milky (IM)	15x5 g
CHF7	ENAMEL PLUS HFO Trial Kit including: 7 syringes 2,5 g with holder 4 dentine: UD2 (A2), UD3 (A3), UD4 (A4), UD5 2 generic enamel: GE2, GE3 1 opalescent enamel: OBN	7x2,5 g

SYRINGE/TIPS HEATER

ENA HEAT Composite Heating Conditioner

It warms up the composite at the temperature of approx. 39°C, making the material more flowable and reducing the stress in the extrusion from the syringe and stratification phases. In this way, this reduces the risk of separation between the matrix and the filler, it improves the biomechanical characteristics of the composite and it avoids bubbles. It is possible to regulate the temperature at 55°C to obtain even more flowable material for cementation. With 6 spaces for syringes and 4 spaces for tips. Also for warming up sodium hypochlorite and anesthetic.



Ref	Description	Q.ty
CHC3	COMPOSITE HEATING CONDITIONER ENA HEAT 110-240V, 50-60 Hz, 500mA Allows to use composite at the ideal temperature of 39°C for modelling and 55°C for cementation	1 pc