

Shiny



Correction of surface tooth geography

Finishing and polishing with rubber point Shiny14, diamond paste Shiny A 3 μ and Shiny B 1 μ and with brushes

Micerium provides everything needed for finishing and polishing of direct esthetic restorations, laboratory veneers, inlays and crowns made in composite.

THE SYSTEM INCLUDES

- Diamond paste: SHINY A 3 μ m, SHINY B 1 μ m
- Aluminium oxide paste SHINY C
- Air Block SHINY G
- Brushes and felts SHINY S, SHINY F
- Pre-polishing rubber point SHINY14 & SHINY33
- Finishing accessories

SURFACE TREATMENT


	Vickers Hardness	Deviat.
Ena HRi microhybrid	55	1,2
HRi microhybrid (a.b.+d.p.)	+10%	0,8
Microfill	18-30	4

Surface treatment: Vickers hardness and its deviation from standard values (a.b.=Air Block • d.p.= Diamond Paste).





Aluminium oxide paste Shiny C with felt wheel; restoration finished and polished under direct light (left) and under polarized light (right)





















Ref	Description	Q.ty
ESKIT	SHINY FINISHING KIT	
	Diamond paste SHINY A 3 micron	2 g
	Diamond paste SHINY B 1 micron	2 g
	Shiny paste SHINY C	35 g
	Air Block SHINY G	7 tips
	SHINY Burs	4 pcs
	Goat hair brushes RA Shiny S/SS	3 pcs
	SHINY 4... brush for posteriors	3 pcs
	SHINY Strips (diamond 5...)	4 pcs
	SHINY Strips (paper M...)	16 pcs
	SHINY Discs	12 pcs
	Felt disc SHINY FD	5 pcs
	Felt wheel SHINY F/SF	5 pcs
	Diamond Rubber point SHINY 14	2 pcs
	Mandrel RA	2 pcs
	Disposable brushes	10 pcs
	Handle for brushes	1 pc
ES+G	POLISHING PASTE SHINY KIT	
	Diamond paste SHINY A 3 micron	2 g
	Diamond paste SHINY B 1 micron	2 g
	Shiny paste SHINY C	35 g
	Goat hair brushes RA Shiny S/SS	3 pcs
	Felt disc SHINY FD & Felt wheel SHINY F/SF	2x3 pcs
	Mandrel RA	1 pc
	Air Block SHINY G	2 tips
	Diamond Rubber point SHINY 14	1 pc
	SHINY 4... brush for posteriors	2 pcs
	SHINY Strips (diamond 5...)	4 pcs
	SHINY Strips (paper M...)	4 pcs
ES-Y+G	POLISHING PASTE SHINY MINI KIT	
	Diamond paste SHINY A 3 micron	0,5 g
	Diamond paste SHINY B 1 micron	0,5 g
	Shiny paste SHINY C	4 g
	Goat hair brushes RA Shiny S	2 pcs
	Mandrel RA	1 pc
	Felt wheel SHINY F	2 pcs
	Diamond Rubber point SHINY 14	1 pc
	Air Block SHINY G	2 tips
	SHINY 4... brush for posteriors	1 pc
	SHINY Strips (diamond 5...)	2 pcs
	SHINY Strips (paper M...)	2 pcs
TSY	FINISHING SHINY KIT FOR LABORATORY	
	Paste SHINY A & Paste SHINY B	2x0,5 g
	Shiny paste SHINY C	4 g
	Tip Glycerine	1 pc
	Paper cone with mandrel	1 pc
	Goat brush Shiny S-HP & Felt Shiny F-HP	2x1 pc
	Carbide Bur FG	1 pc
	Green carborundum bur HP	1 pc
	Diamond burs HP (D717-D725-D727)	3 pcs
	Diamond disc 2 sides 0,15mm D751	1 pc

ENA SHINY REFILL

Ref	Description	Q.ty
ESG7	SHINY G glycerine Air Block in tips 0,5 ml	7 pcs
ESG	SHINY G glycerine Air Block in syringe	2 ml
ESG10	SHINY G glycerine Air Block in syringe	10 ml
ESA	SHINY A 3 micron Diamond paste for pre-polishing to be used with Shiny S	2 g
ESB	SHINY B 1 micron Diamond paste for polishing to be used with Shiny S	2 g
ESC	SHINY C aluminium oxide paste for shiny to be used with Shiny F	35 g
ESS	SHINY S goat brush wheel RA	12 pcs
ESSHP	SHINY S goat brush wheel HP	12 pcs
CSS	SHINY SS goat brush wheel RA Mini for posterior	12 pcs
ESF	SHINY F felt wheel 2 mm to be used with ES6	30 pcs
ESFD	SHINY FD interproximal felt disc 0,5 mm	100 pcs
CSF	SHINY SF cotton felt wheel RA (17x10 mm)	12 pcs
ESFHP	SHINY FHP felt wheel	12 pcs
ES13	SHINY rubber points for pre-polishing	12 pcs
ES33	SHINY brown rubber points	12 pcs
CS014	SHINY14 Diamond rubber points for pre-polishing composite	6 pcs
CS015	SHINY15 Diamond rubber for polishing composite	3 pcs
CS041	SHINY41 Impregnated Brush for posteriors "Point"	6 pcs
CS042	SHINY42 Impregnated Brush for posteriors "Mini Cup"	6 pcs
CS043	SHINY43 Impregnated Brush for posteriors "Cup"	6 pcs





Ref	Description	Q.ty	
COSSHINYCON	Paper cones HP	25 pcs	
COSSHINYCONM	Mandrel for paper cones Shiny	1 pc	
ES6	SHINY Mandrel for felt wheels	6 pcs	
ESFM	SHINY Mandrel for discs	1 pc	
ESDF	SHINY DISCS extra fine brown grain	100 pcs	
ESDX	SHINY DISCS fine blue grain	100 pcs	
ESDM	SHINY DISCS medium green grain	100 pcs	
ESDC	SHINY DISCS coarse ivory grain	100 pcs	
ESSF	SHINY STRIPS fine (blue) extrafine (brown) grain	100 pcs	
ESSFS	SHINY STRIPS fine-extra fine narrow grain	100 pcs	
ESSM	SHINY STRIPS medium (green) coarse (ivory) grain	100 pcs	
ESSMS	SHINY STRIPS medium-coarse narrow grain	100 pcs	
COSSHINY1	Handle for brushes	1 pc	
COSSHINY50	Brush refill	50 pcs	
ES25	Tips ø 0,9 mm for Glycerine & Flow	25 pcs	
CS050-2R	Diamond strips 120 mm x 2 mm fine red	5 pcs	
CS050-2Y	Diamond strips 120 mm x 2 mm extra fine yellow	5 pcs	
CS050-4R	Diamond strips 120 mm x 4 mm fine red	5 pcs	
CS050-4Y	Diamond strips 120 mm x 4 mm extra fine yellow	5 pcs	